

# Respirator mask

FFP2 & FFP3



## FFP2, FFP2 Color and FFP 3

FFP2 masks protect against harmful water- and oil-based particles. FFP3 masks also additionally protect against carcinogenic substances, radioactive particles, biological agents of risk group 3 and enzymes.

Our respirators are also characterized by particularly pleasant and effective material and high wearing comfort.

### Mouldable nose bridge

Tight closure is ensured by the moldable nose clip, which is incorporated into the mask. The inside of the nose clip is reinforced with comfortable foam (FFP3 only).

### Made in EU

Both mask types are manufactured in the Czech Republic (EU). This guarantees long-term good and sustainable quality, as well as short delivery routes.

### Light breathing

The masks are characterized by their ultra-thin soft filter material (made in CZ).

### Closure at back of head

An enclosed S-shaped closure cap allows the mask to be fixed to the back of the head.

### Comfortable loops


The soft masks have comfortable retaining straps and allows comfortable wearing.

### Individually packed

In each box there are 5 FFP3 or 10 FFP2 masks, individually packed.

### Certified mask

The white FFP2 and FFP3 masks are certified according to EN 149:2001 + A1:2009 (PPE) and EN 14683:2019 + AC:2019 (OR approval).



**The FFP2 mask fulfills all requirements to be worn in the operating room. The mask protects the wearer and the patient at the same time.**



## Certified according to EN 149:2001 + A1:2009 and EN 14683:2019 + AC:2019

Our HERZmed FFP2 and FFP3 mask in white is certified not only for self-protection according to EN 14683:2019+AC:2019 (PPE), but also as a medical device according to EN 149:2001+A1:2009 (external protection).

**NOW MEDICAL  
CERTIFICATED AS  
OP-MASK**

### Comfort that convinces!

Our masks are made of high quality filter material, which is particularly comfortable to wear. The soft and stretchy ear loops allow wearing without unpleasant pressure points even for several hours.

- |                            |                                |
|----------------------------|--------------------------------|
| ✓ <b>Comfortable loops</b> | ✓ <b>With closing flap</b>     |
| ✓ <b>No pressure marks</b> | ✓ <b>Mouldable nose bridge</b> |
| ✓ <b>Soft material</b>     | ✓ <b>Light breathing</b>       |



### 100 % Made in EU

A very important issue in the production of face, mouth and respiratory masks is undoubtedly the material used in the product. Therefore, our Czech manufacturing company uses only filter material made in the Czech Republic.



- |                                 |                                       |
|---------------------------------|---------------------------------------|
| ✓ <b>Short supply chains</b>    | ✓ <b>No unnecessary fleece layers</b> |
| ✓ <b>CO<sub>2</sub> Savings</b> | ✓ <b>High quality</b>                 |
| ✓ <b>No sea containers</b>      | ✓ <b>Fair prices</b>                  |

### Highlights

- Certified masks
- 100% Made in EU
- High quality material
- Comfortable loops
- Individually packed
- Moldable nose clip
- Easy breathing
- Individual logo printing possible (on request)





## Technical specifications

### FFP3 mask (Box with 5 pieces):

**Size box:** 130 x 65 x 197 mm (L x B x H)  
**Size mask standard:** 145 x 90 mm (L x B)  
**Weight Karton:** 80 g  
**Coloring:** Weiß

### FFP2 mask (Box with 10 pieces):

**Size box:** 130 x 65 x 197 mm (L x B x H)  
**Size mask standard:** 155 x 105 mm (L x B)  
**Size mask L:** 155 x 95 mm (L x B)  
**Size mask M:** 145 x 90 mm (L x B)  
**Weight box:** 100 g  
**Coloring:** White/Color

### Material:


Inner layer: Spunbonded nonwoven (PP)  
 Middle filter layer: Meltblown fleece  
 Outer layer: Spunbond (PP)

### Sterility:









No

## Order information

FFP3

Art.-No.	Description	Image
M6320002	FFP3 mask, size „standard“, 5 pieces, white	

FFP2

Art.-No.	Description	Image
311.10011	FFP2 mask, size „standard“, 10 pieces, white	
311.10022	FFP2 mask, size „M“, 10 pieces, white	
311.10023	FFP2 mask, size „L“, 10 pieces, white	
311.10014	FFP2 mask, size „standard“, 10 pieces, grey	
311.10015	FFP2 mask, size „standard“, 10 pieces, light blue	
311.10016	FFP2 mask, size „standard“, 10 pieces, ochre	
311.10017	FFP2 mask, size „standard“, 10 pieces, pink	
311.10018	FFP2 mask, size „standard“, 10 pieces, light pink	
311.10019	FFP2 mask, size „standard“, 10 pieces, red	
311.10020	FFP2 mask, size „standard“, 10 pieces, black	
311.10021	FFP2 mask, size „standard“, 10 pieces, turquoise	

● With TYPE IIR certification

## Order and advice by:

