

# Emergency containers



## The material makes the difference *WaterStop® - for demanding applications*

We are pleased to introduce our new product line WaterStop. We have been supplying a very successful range of emergency containers for years. Our primary aim is not to invent new models, but to improve our existing range through innovation, higher quality high-tech materials and manufacturing processes. Our aim is to respond to existing requirements and needs of professional users.

**WaterStop products are designed for professional use.**



**Water repellent zips**

The particularly robust and special zips are water-repellent and therefore also suitable for use in the rain and allow the storage of dressing material that is susceptible to moisture.


- Water-repellent material
- Easier cleaning of conventional soiling
- Improved durability
- Easy opening and closing



**PU-coated material**

The polyester fabric is coated twice with a high-quality PU coating. This has numerous advantages over the conventional PVC coating.

- Significantly lower weight
- Higher load capacity in daily use
- Better temperature resistance, especially at low the material remains flexible
- Improved durability
- More environmentally friendly than conventional PVC material



**WaterStop® Technology with TEFLON®**

Polyester fabric with TEFLON® coating, which refines ordinary polyester into a unique material. Water and oil can no longer penetrate the fabric and are repelled on the surface.

- Water repellent material
- Easier to clean from conventional dirt
- Significantly reduced weight
- Improved durability



## Attention to detail

Emergency containers face tough challenges in the field. For this reason, HERZmed has set itself the task of making it its business to support the user with high-quality and robust materials. With the new WaterStop series, an even more resistant material is now available. In order to provide you with quick and safe access to the products in an emergency, you also have the following benefits:

### WaterStop Emergency Backpack PRO

#### Sturdy grab handle

The ergonomically shaped carrying handle makes it comfortable to carry even heavy emergency containers.

#### Extra strong zippers

In extreme situations, every move must be perfect. Our innovative HERZmed® zippers are particularly stable and provide optimum guidance. Thanks to the rubber coating and the relief in the shape of the HERZmed heart, the zippers fit perfectly in the hand and thus ensure quick access to the contents.

#### Water repellent zip

The use of water-repellent zips on all exposed areas provides optimal conditions for storing dressing material that is susceptible to moisture. High quality ensures long life and easy opening and closing.

#### PVC viewing tray

The transparent viewing compartment at the front is protected by a water-repellent zip in the rain and is also suitable for sensitive documents.

#### Quality down to the smallest detail

A highlight of our WaterStop line are the elaborately sewn-in reflective strips. Compared to conventional reflective strips sewn on openly, these are much better protected and therefore more durable. We have also replaced all buckles and stressed eyelets with more robust versions.

#### Own logo possible

You have the option of attaching your own logo, which can be quickly changed using Velcro, to each backpack.

#### Customisable

A printable and changeable reflective strip on the front offers you numerous possibilities to fulfil your individual wishes.



## Attention to detail

Emergency containers face tough challenges in the field. For this reason, HERZmed has set itself the task of making it its business to support the user with high-quality and robust materials. With the new WaterStop series, an even more resistant material is now available. In order to provide you with quick and safe access to the products in an emergency, you also have the following benefits:

### **PREMIUM Emergency Backpack X1 YELLOW PLANE**

#### **Sturdy grab handle**

The ergonomically shaped carrying handle makes it comfortable to carry even heavy emergency containers.

#### **Extra strong zippers**

In extreme situations, every move must be perfect. Our innovative HERZmed® zippers are particularly stable and provide optimum guidance. Thanks to the rubber coating and the relief in the shape of the HERZmed heart, the zippers fit perfectly in the hand and thus ensure quick access to the contents.

#### **Water repellent zip**

The use of water-repellent zips on all exposed areas provides optimal conditions for storing dressing material that is susceptible to moisture. High quality ensures long life and easy opening and closing.

#### **PVC viewing tray**

The transparent viewing compartment at the front is protected by a water-repellent zip in the rain and is also suitable for sensitive documents.

#### **O2 bottle fixation**

With practical interior partitioning, extra compartment for O2 bottles up to 2 litres. The fixation velcro is removable.

#### **Own logo possible**

You have the option of attaching your own logo, which can be quickly changed using Velcro, to each backpack.

#### **Edge protection**

The robust edge protection on the corners ensures that the emergency backpack lasts longer.

#### **Reflective stripes**

Two sewn-in reflective strips provide more safety and visibility during use.

#### **Detachable shoulder straps**

The shoulder straps of the Premium Emergency Backpack X1 can be temporarily removed and then reattached as needed.





## Attention to detail

Emergency containers face tough challenges in the field. For this reason, HERZmed has set itself the task of making it its business to support the user with high-quality and robust materials. With the new WaterStop series, an even more resistant material is now available. In order to provide you with quick and safe access to the products in an emergency, you also have the following benefits:

### WaterStop Emergency Bag ULTRA

#### Water repellent zip

The use of water-repellent zips on all exposed areas provides optimal conditions for storing dressing material that is susceptible to moisture. High quality ensures long life and easy opening and closing.

#### Metal eyelets for shoulder strap

There are two sturdy metal eyelets on the left and right. The enclosed shoulder strap can be attached to these.

#### Side compartments

The WaterStop Emergency Bag Ultra has two practical side pockets.

#### Front compartment

The emergency bag scores with its many compartments and various divisions for maximum order. In addition to the two side compartments, it also has a large front pocket.

#### Reflective stripes

Two sewn-in reflective strips provide more safety and visibility during use.

#### Velcro fastener as grab handle

A flexible handle with Velcro fastener makes the WaterStop Emergency Bag Ultra comfortable to carry.

#### Inner compartments

The bag also has two inner compartments. The lower of the two is equipped with 7 elastic holding loops (e.g. for Guedel tubes).

#### Fixed distribution

The main compartment is characterised by a fixed inner division.

#### Own logo possible

You have the option of attaching your own logo, which can be quickly changed using Velcro, to each backpack.

#### Customisable

A printable and changeable reflective strip on the front offers you numerous possibilities to fulfil your individual wishes.



## WaterStop Emergency Backpack PROFI

Measures: 65 × 42 × 23 cm

### Highlights

- Teflon® SHIELD+ material with beading effect
- Water repellent, dirt repellent
- Special, robust zips for use in the rain
- Grey lining for good visibility of contents
- Sewn-in reflective stripes with edge protection
- Practical outer pockets
- 4 module pockets included
- O2 bottle attachment

Art.-No.	Article description	Material	Colour
N180012	WaterStop Emergency Backpack PROFI	Plane	●
N180009	WaterStop Emergency Backpack PROFI	Nylon	●
N180066	WaterStop Emergency Backpack PROFI	Plane	●
N180068	WaterStop Emergency Backpack PROFI	Nylon	●
N180067	WaterStop Emergency Backpack PROFI	Plane	●
N180069	WaterStop Emergency Backpack PROFI	Nylon	●

Empty weight nylon-version: approx. 3.5 kg, plane-version approx. 5.0 kg



## WaterStop Emergency Backpack PRO

Measures: 46 × 25 × 20 cm

### Highlights

- Teflon® SHIELD+ material with beading effect
- Water repellent, dirt repellent
- Special, robust zips for use in the rain
- Grey lining for good visibility of contents
- Sewn-in reflective stripes
- 2 side and front pockets (lower pocket suitable for AED)
- Including 5 module pockets

Art.-No.	Article description	Material	Colour
N180005	WaterStop Emergency Backpack PRO	Plane	●
N180002	WaterStop Emergency Backpack PRO	Nylon	●
N180006	WaterStop Emergency Backpack PRO	Plane	●
N180003	WaterStop Emergency Backpack PRO	Nylon	●
N180004	WaterStop Emergency Backpack PRO	Nylon	●
N180070	WaterStop Emergency Backpack PRO	Plane	●
N180072	WaterStop Emergency Backpack PRO	Nylon	●
N180071	WaterStop Emergency Backpack PRO	Plane	●
N180073	WaterStop Emergency Backpack PRO	Nylon	●

Empty weight nylon-version: approx. 1.5 kg, plane-version approx. 2.7 kg





## WaterStop Emergency Backpack SMALL RED

Measures: 43 × 35 × 15 cm

### Highlights

- Teflon® SHIELD+ material with beading effect
- Water repellent, dirt repellent
- Special, robust zips for use in the rain
- Grey lining for good visibility of contents
- Sewn-in reflective stripes
- Front pocket
- Includes 3 module pockets

Art.-No.	Article description	Material	Colour
N180001	WaterStop Emergency Backpack SMALL RED	Nylon	●

Empty weight nylon-version: approx. 0.9 kg



## WaterStop Emergency Bag FREE

Measures: 44 × 30 × 18 cm

### Highlights

- Teflon® SHIELD+ material with beading effect,
- Water-repellent, dirt-repellent
- Special, robust zips for use in the rain
- Grey lining for good visibility of contents
- Sewn-in reflective stripes
- Structure: Main compartment with padded, 3-compartment inner divider (removable), front compartment, 2 side compartments

Art.-No.	Article description	Material	Colour
N180014	WaterStop Emergency Bag FREE	Nylon	●
N180015	WaterStop Emergency Bag FREE	Nylon	●
N180016	WaterStop Emergency Bag FREE	Nylon	●
N180080	WaterStop Emergency Bag FREE	Nylon	●
N180081	WaterStop Emergency Bag FREE	Nylon	●

Empty weight nylon-version: approx. 0.88 kg



## WaterStop Emergency Bag MULTI

Measures: 45 × 24 × 20 cm

### Highlights

- Teflon® SHIELD+ material with beading effect
- Water repellent, dirt repellent
- Special, robust zips for use in the rain
- Grey lining for good visibility of contents
- Sewn-in reflective stripes
- Structure: Main compartment with fixed inner divider, 2 side compartments

Art.-No.	Article description	Material	Colour
N180017	WaterStop Emergency Bag MULTI	Nylon	Red
N180018	WaterStop Emergency Bag MULTI	Nylon	Blue
N180019	WaterStop Emergency Bag MULTI	Nylon	Black
N180083	WaterStop Emergency Bag MULTI	Nylon	Yellow

Empty weight nylon-version: approx. 0.97 kg



## WaterStop Emergency Bag ULTRA


Measures: 44 × 34 × 30 cm

### Highlights

- Teflon® SHIELD+ material with beading effect
- Water repellent, dirt repellent
- Special, robust zips for use in the rain
- Grey lining for good visibility of contents
- Sewn-in reflective stripes
- Structure: Main compartment with fixed inner divider, inner compartment with 7 elastic holders (e.g. for Guedel tubes) Front compartment, 2 side compartments

Art.-No.	Article description	Material	Colour
N180023	WaterStop Emergency Backpack ULTRA	Plane	Red
N180020	WaterStop Emergency Backpack ULTRA	Nylon	Red
N180024	WaterStop Emergency Backpack ULTRA	Plane	Blue
N180021	WaterStop Emergency Backpack ULTRA	Nylon	Blue
N180022	WaterStop Emergency Backpack ULTRA	Nylon	Black
N180084	WaterStop Emergency Backpack ULTRA	Nylon	Green
N180085	WaterStop Emergency Backpack ULTRA	Nylon	Yellow

Empty weight nylon-version: approx. 1.7 kg





## WaterStop Ampoule Kit PRO

Measures: 29 × 24 × 8,5 cm

### Highlights

- Teflon® SHIELD+ material with beading effect
- Water repellent, dirt repellent
- Robust zips
- Grey lining for good visibility of contents
- Sewn-in reflective stripes

Art.-No.	Article description	Material	Colour
N180031	WaterStop Ampoule Kit PRO	Plane	●
N180034	WaterStop Ampoule Kit PRO	Nylon	●
N180075	WaterStop Ampoule Kit PRO	Plane	●

Empty weight nylon-version: approx. 0,5 kg, plane-version approx. 0.6 kg



## WaterStop Ampoule Kit LIGHT

Measures: 13 × 18 × 6 cm

### Highlights

- Teflon® SHIELD+ material with beading effect,
- Water repellent, dirt repellent
- Robust zips
- Grey inner lining for good visibility of contents

Art.-No.	Article description	Material	Colour
N180028	WaterStop Ampoule Kit LIGHT	Plane	●
N180025	WaterStop Ampoule Kit LIGHT	Nylon	●
N180026	WaterStop Ampoule Kit LIGHT	Nylon	●
N180030	WaterStop Ampoule Kit LIGHT	Plane	●
N180077	WaterStop Ampoule Kit LIGHT	Plane	●
N180079	WaterStop Ampoule Kit LIGHT	Nylon	●

Empty weight nylon-version: approx. 0.2 kg, plane-version approx. 0.28 kg



## WaterStop O2-Carrier Bag PROFI BLUE

Measures: 49 x 24 x 17 cm

### Highlights

- Suitable for O2 bottles up to 2 litres
- Storage space through additional side pockets (1x outside; 2x inside)
- Holder for carrying and shoulder strap (detachable)
- Carrying strap with metal clips
- High-quality metal quick-release fastener on the back
- Hinged clear window
- Sewn-in white reflective stripes
- Fits perfectly in the emergency vehicle
- Can be attached to an emergency stretcher

Art.-No.	Article description	Material	Colour
N190308	WaterStop O2-Carrier Bag PROFI BLUE	Plane	Blue

Empty weight planen-version: approx. 1.57 kg



## O2-Carrier Bag BLACK

Measures: 48 x 22 x 11 cm

### Highlights

- Durable material
- Water repellent
- Comfortable carrying strap
- O2 bottle up to 2 litres fixable

Art.-No.	Article description	Material	Colour
M660013	O2-Carrier Bag BLACK	Nylon	Black

Empty weight nylon-version: approx. 0.87





## O2-Carrier Bag BLUE PLANE

Measures: 48 × 22 × 11 cm

### Highlights

- Suitable for O2 cylinders up to 2 litres in combination with a pressure reducer
- Storage space due to additional side pocket
- One additional module pocket inside
- Holder for carrying and shoulder strap (detachable)
- Hinged clear-view window allows a direct view of the pressure reducer and also quick access

Art.-No.	Article description	Material	Colour
N190306	O2-Carrier Bag BLUE PLANE	Plane	<span style="color: blue;">●</span>

Empty weight nylon-version: approx. 1.4 kg



## Emergency Bag FLEX PLANE

Measures: 28 × 32 × 18 cm

### Highlights

- Compartmentation: 1 main compartment with fold-out interior, three compartments that can be closed with a zip, many elastic straps (for e.g. Guedel tubes) and a modular pocket (ideal for mobile dressing material)
- Tear-resistant tarpaulin material, water-repellent
- Grey inner lining for good visibility of contents
- Optionally with hooks for wall mounting
- For stationary and ambulant use

Art.-No.	Article description	Material	Colour
N180035	WaterStop Emergency Bag FLEX PLANE	Plane	<span style="color: red;">●</span>
N180064	WaterStop Emergency Bag FLEX PLANE	Plane	<span style="color: green;">●</span>
N180065	WaterStop Emergency Bag FLEX PLANE	Plane	<span style="color: yellow;">●</span>

Empty weight nylon-version: approx. 1.55 kg



## Wallhooks

### Art.-Nr.

N180044

- For wall mounting of the emergency bag FLEX



## Rescue sheet PRO PLANE RED SAFE with footmuff

Measures: 200 × 70 cm

### Highlights

- Footmuff for stable and safe transport
- Can be used as a rescue seat
- 10 stable handles
- Load capacity of up to 150 kg
- DEKRA Germany tested (complies with DIN EN 1865)

Art.-No.	Article description	Material	Colour
N190312	Rescue sheet PRO PLANE RED SAFE	Plane	●



## Logbook A5 Organiser Stapler PLANE

Measures: 28 × 22 × 5 cm

### Highlights

- Tried and tested tarpaulin material
- Flexible division
- Suitable for DIN A5 documents
- With practical stapler and key ring

Art.-No.	Article description	Material	Colour
N180043	Logbook A5 Organiser Stapler PLANE	Plane	●
N180042	Logbook A5 Organiser Stapler PLANE	Plane	●

Empty weight: approx. 0.43 kg





## Logbook A5 Organiser PLANE

Measures: 28 × 22 × 5 cm

### Highlights

- Tried and tested tarpaulin material
- Flexible division
- Suitable for DIN A5 documents
- With key ring

Art.-No.	Article description	Material	Colour
N180041	Logbook A5 Organiser PLANE	Plane	Red
N180040	Logbook A5 Organiser PLANE	Plane	Blue

Empty weight: approx. 0.43 kg



## MULTI Operation Organiser A4 PLANE

Measures: 28 × 37 × 6 cm

### Highlights

- Eye-catching white reflective stripes
- Tried and tested tarpaulin material
- Flexible division
- Triage colour division
- Suitable for DIN A4 documents
- With practical stapler and key ring

Art.-No.	Article description	Material	Colour
N180038	MULTI Operation Organiser A4 PLANE	Plane	Red
N180037	MULTI Operation Organiser A4 PLANE	Plane	Blue

Empty weight: approx. 0.85 kg



## PREMIUM Emergency Backpack X1 YELLOW PLANE

Measures: 52 × 30 × 25 cm

### Highlights

- Tried and tested tarpaulin material
- Emergency bag/backpack combination
- Includes 7 colour-coded modular bags with viewing window
- Sewn-in reflective strips with edge protection
- Removable O2 bottle attachment

Art.-No.	Article description	Material	Colour
N190189-V2	PREMIUM Emergency Backpack X1 YELLOW PLANE	Plane	Red
N190171-V2	PREMIUM Emergency Backpack X1 YELLOW PLANE	Plane	Blue
N190441	PREMIUM Emergency Backpack X1 YELLOW PLANE	Plane	Green
N190442	PREMIUM Emergency Backpack X1 YELLOW PLANE	Plane	Yellow

Empty weight planen-version: approx. 6.98 kg



## Emergency Backpack PROFi

Measures: 28 × 5 × 22 cm

### Highlights

- Made from durable material
- Practical outer pockets
- 4 module pockets with viewing window
- Yellow reflective stripes
- O2 bottle attachment for up to 2 litres

Art.-No.	Article description	Material	Colour
N190115	Emergency Backpack PROFi	Plane	Red
N190117	Emergency Backpack PROFi	Plane	Blue
N190114	Emergency Backpack PROFi	Nylon	Red

Empty weight nylon-version: approx. 3.79 kg, plane-version approx. 5.15 kg





## Emergency Backpack YELLOW LARGE

Measures: 36 × 47 × 26 cm

### Highlights

- Made of durable material - polyester nylon or tarpaulin
- Practical outer pockets
- 5 colour coded module pockets with viewing window
- Yellow reflective stripes
- Freely divisible interior
- Clear view pocket

Art.-No.	Article description	Material	Colour
N190134	Emergency Backpack YELLOW LARGE	Plane	●
N190132	Emergency Backpack YELLOW LARGE	Nylon	●
N190244	Emergency Backpack YELLOW LARGE	Plane	●
N190243	Emergency Backpack YELLOW LARGE	Nylon	●
N190133	Emergency Backpack YELLOW LARGE	Nylon	●
N190445	Emergency Backpack YELLOW LARGE	Nylon	●
N190446	Emergency Backpack YELLOW LARGE	Nylon	●

Empty weight nylon-version: approx. 1.5 kg, plane-version approx. 2.5 kg



## Emergency Bag YELLOW SMALL

Measures: 37 × 21 × 27 cm

### Highlights

- Made from durable material
- Practical outer pockets
- Three-part interior

Art.-No.	Article description	Material	Colour
N190140	Emergency Bag YELLOW SMALL	Nylon	●
N190203	Emergency Bag YELLOW SMALL	Nylon	●
N190447	Emergency Bag YELLOW SMALL	Nylon	●
N190448	Emergency Bag YELLOW SMALL	Nylon	●

Empty weight nylon-version: approx. 0.77



## Standard Module Bags



Art.-No.	Article description	Measures	Material	Colour
N190124	Module Bag S Blue	24 x 12 x 9 cm	Plane	Blue
N190127	Module Bag S Blue	24 x 12 x 9 cm	Nylon	Blue
N190128	Module Bag S Red	24 x 12 x 9 cm	Nylon	Red
N190129	Module Bag S Lemon	24 x 12 x 9 cm	Nylon	Yellow
N190130	Module Bag S Orange	24 x 12 x 9 cm	Nylon	Orange
N190131	Module Bag S Green	24 x 12 x 9 cm	Nylon	Green
N190087	Module Bag S SET	24 x 12 x 9 cm	Nylon	1 x each



Art.-No.	Article description	Measures	Material	Colour
N190123	Module Bag M Red	32 x 12 x 5 cm	Plane	Red
N190125	Module Bag M Green	32 x 12 x 5 cm	Nylon	Green



Art.-No.	Article description	Measures	Material	Colour
N190126	Module Bag L Blue	23 x 20 x 10 cm	Nylon	Blue
N190122	Module Bag L Green	23 x 20 x 10 cm	Plane	Green
N190291	Module Bag L Black	23 x 20 x 10 cm	Plane	Black

## For example suitable for





## WaterStop Module Bags




Art.-No.	Article description	Measures	Material	Colour
N180086	WaterStop Module Bag S Blue	24 x 12 x 9 cm	Nylon	Blue
N180087	WaterStop Module Bag S Green	24 x 12 x 9 cm	Nylon	Green
N180088	WaterStop Module Bag S Red	24 x 12 x 9 cm	Nylon	Red
N180089	WaterStop Module Bag S Lemon	24 x 12 x 9 cm	Nylon	Yellow
N180090	WaterStop Module Bag S Black	24 x 12 x 9 cm	Nylon	Black



Art.-No.	Article description	Measures	Material	Colour
N180107	WaterStop Module Bag S Blue	24 x 12 x 9 cm	Plane	Blue
N180108	WaterStop Module Bag S Green	24 x 12 x 9 cm	Plane	Green
N180109	WaterStop Module Bag S Red	24 x 12 x 9 cm	Plane	Red
N180110	WaterStop Module Bag S Lemon	24 x 12 x 9 cm	Plane	Yellow
N180111	WaterStop Module Bag S Black	24 x 12 x 9 cm	Plane	Black



Art.-No.	Article description	Measures	Material	Colour
N180093	WaterStop Module Bag M Blue	32 x 12 x 5 cm	Nylon	Blue
N180094	WaterStop Module Bag M Green	32 x 12 x 5 cm	Nylon	Green
N180095	WaterStop Module Bag M Red	32 x 12 x 5 cm	Nylon	Red
N180096	WaterStop Module Bag M Lemon	32 x 12 x 5 cm	Nylon	Yellow
N180097	WaterStop Module Bag M Black	32 x 12 x 5 cm	Nylon	Black



Art.-No.	Article description	Measures	Material	Colour
N180100	WaterStop Module Bag L Blue	23 x 20 x 10 cm	Nylon	Blue
N180101	WaterStop Module Bag L Green	23 x 20 x 10 cm	Nylon	Green
N180102	WaterStop Module Bag L Red	23 x 20 x 10 cm	Nylon	Red
N180103	WaterStop Module Bag L Lemon	23 x 20 x 10 cm	Nylon	Yellow
N180104	WaterStop Module Bag L Black	23 x 20 x 10 cm	Nylon	Black

## For example suitable for



HERZmed WaterStop  
Emergency Backpack PROFI



HERZmed WaterStop  
Emergency Backpack PRO



HERZmed Emergency  
Backpack YELLOW Large



HERZmed WaterStop  
Emergency Backpack SMALL



HERZmed WaterStop  
Emergency Bag FLEX

## Module Bags PROFI



Art.-No.	Article description	Measures	Material	Colour
N190370	WaterStop Module Bag PROFI XS Blue	16 x 15 x 6 cm	Plane	Blue
N190371	WaterStop Module Bag PROFI XS Green	16 x 15 x 6 cm	Plane	Green
N190372	WaterStop Module Bag PROFI XS Red	16 x 15 x 6 cm	Plane	Red
N190373	WaterStop Module Bag PROFI XS Lemon	16 x 15 x 6 cm	Plane	Yellow
N190374	WaterStop Module Bag PROFI XS Black	16 x 15 x 6 cm	Plane	Black
N190375	WaterStop Module Bag PROFI XS Grey	16 x 15 x 6 cm	Plane	Grey



Art.-No.	Article description	Measures	Material	Colour
N190377	WaterStop Module Bag PROFI S Blue	22 x 15 x 7 cm	Plane	Blue
N190378	WaterStop Module Bag PROFI S Green	22 x 15 x 7 cm	Plane	Green
N190379	WaterStop Module Bag PROFI S Red	22 x 15 x 7 cm	Plane	Red
N190380	WaterStop Module Bag PROFI S Lemon	22 x 15 x 7 cm	Plane	Yellow
N190381	WaterStop Module Bag PROFI S Black	22 x 15 x 7 cm	Plane	Black
N190382	WaterStop Module Bag PROFI S Grey	22 x 15 x 7 cm	Plane	Grey



Art.-No.	Article description	Measures	Material	Colour
N190384	WaterStop Module Bag PROFI M Blue	25 x 12 x 6 cm	Plane	Blue
N190385	WaterStop Module Bag PROFI M Green	25 x 12 x 6 cm	Plane	Green
N190386	WaterStop Module Bag PROFI M Red	25 x 12 x 6 cm	Plane	Red
N190387	WaterStop Module Bag PROFI M Lemon	25 x 12 x 6 cm	Plane	Yellow
N190388	WaterStop Module Bag PROFI M Black	25 x 12 x 6 cm	Plane	Black
N190389	WaterStop Module Bag PROFI M Grey	25 x 12 x 6 cm	Plane	Grey



Art.-No.	Article description	Measures	Material	Colour
N190391	WaterStop Module Bag PROFI L Blue	37 x 14 x 10 cm	Plane	Blue
N190392	WaterStop Module Bag PROFI L Green	37 x 14 x 10 cm	Plane	Green
N190393	WaterStop Module Bag PROFI L Red	37 x 14 x 10 cm	Plane	Red
N190394	WaterStop Module Bag PROFI L Lemon	37 x 14 x 10 cm	Plane	Yellow
N190395	WaterStop Module Bag PROFI L Black	37 x 14 x 10 cm	Plane	Black
N190396	WaterStop Module Bag PROFI L Grey	37 x 14 x 10 cm	Plane	Grey

## For example suitable for



HERZmed WaterStop  
Emergency Backpack PROFI



HERZmed WaterStop  
Emergency Backpack PRO



HERZmed Emergency  
Backpack YELLOW Large



## Module Bags X1



Art.-No.	Article description	Measures	Material	Colour
N190349	WaterStop Module Bag X1 S Blue	21 x 15 x 10 cm	Plane	Blue
N190350	WaterStop Module Bag X1 S Green	21 x 15 x 10 cm	Plane	Green
N190351	WaterStop Module Bag X1 S Red	21 x 15 x 10 cm	Plane	Red
N190352	WaterStop Module Bag X1 S Lemon	21 x 15 x 10 cm	Plane	Yellow
N190354	WaterStop Module Bag X1 S Grey	21 x 15 x 10 cm	Plane	Grey



Art.-No.	Article description	Measures	Material	Colour
N190342	WaterStop Module Bag X1 L Blue	30 x 13 x 10 cm	Plane	Blue
N190343	WaterStop Module Bag X1 L Green	30 x 13 x 10 cm	Plane	Green
N190344	WaterStop Module Bag X1 L Red	30 x 13 x 10 cm	Plane	Red
N190345	WaterStop Module Bag X1 L Lemon	30 x 13 x 10 cm	Plane	Yellow
N190347	WaterStop Module Bag X1 L Grey	30 x 13 x 10 cm	Plane	Grey



Module Bag L in X1 Emergency Backpack



Module Bag S in X1 Emergency Backpack

## Suitable for



HERZmed PREMIUM  
Emergency Backpack X1

## Ordering information

Image	Unit	PU	Art. No.	Description	Volume	Measures	SRP
	Piece	2	N180012 (red) N180066 (green on request) N180067 (yellow)	WaterStop Emergency Backpack PROFI (PLANE)	approx. 67,1 L	approx. 65 x 42 x 23 cm	180,00 €**
	Piece	2	N180009 (red) N180068 (green on request) N180069 (yellow)	WaterStop Emergency Backpack PROFI (NYLON)	approx. 67,1 L	approx. 65 x 42 x 23 cm	135,00 €**
	Piece	4	N180005 (red) N180006 (Blue) N180070 (green on request) N180071 (yellow)	WaterStop Emergency Backpack PRO (PLANE)	approx. 30 L	approx. 46 x 25 x 20 cm	85,00 €**
	Piece	10	N180002 (red) N180003 (Blue) N180004 (Black) N180072 (green on request) N180073 (yellow)	WaterStop Emergency Backpack PRO (NYLON)	approx. 30 L	approx. 46 x 25 x 20 cm	65,00 €**
	Piece	20	N180001 (red)	WaterStop Emergency Backpack SMALL (NYLON)	approx. 22,8 L	approx. 43 x 35 x 15 cm	40,00 €**
	Piece	12	N180014 (red) N180015 (red) N180016 (Black) N180080 (green on request) N180081 (yellow)	WaterStop Emergency Bag FREE (NYLON)	approx. 16,6 L	approx. 44 x 30 x 18 cm	45,00 €**
	Piece	25	N180017 (red) N180018 (Blue) N180019 (Black) N180083 (yellow)	WaterStop Emergency Bag MULTI (NYLON)	approx. 18,5 L	approx. 45 x 24 x 20 cm	45,00 €**
	Piece	6	N180023 (red) N180024 (Blue)	WaterStop Emergency Bag ULTRA (PLANE)	approx. 36,4 L	approx. 44 x 34 x 30 cm	85,00 €**
	Piece	10	N180020 (red) N180021 (Blue) N180022 (Black) N180084 (green on request) N180085 (yellow)	WaterStop Emergency Bag ULTRA (NYLON)	approx. 36,4 L	approx. 44 x 34 x 30 cm	65,00 €**
	Piece	20	N180031 (red) N180075 (yellow)	WaterStop Ampoule Kit PRO (PLANE)	approx. 5,9 L	approx. 29 x 24 x 8,5 cm	45,00 €**
	Piece	20	N180034 (red)	WaterStop Ampoule Kit PRO (NYLON)	approx. 5,9 L	approx. 29 x 24 x 8,5 cm	35,00 €**
	Piece	50	N180028 (red) N180030 (Black) N180077 (yellow)	WaterStop Ampoule Kit LIGHT (PLANE)	approx. 1,4 L	approx. 13 x 18 x 6 cm	15,00 €**
	Piece	50	N180025 (red) N180026 (Blue) N180079 (yellow)	WaterStop Ampoule Kit LIGHT (NYLON)	approx. 1,4 L	approx. 13 x 18 x 6 cm	12,50 €**
	Piece	10	M660013 (Black)	O2 - Carrier bag MOBIL (NYLON)	approx. 11 L	approx. 48 x 22 x 11 cm	40,00 €**
	Piece	4	N190306 (Blue)	Oxygen Carrying Bag Model 2020 (PLANE)	approx. 11,6 L	approx. 48 x 22 x 11 cm	70,00 €**
	Piece	6	N190308 (Blue)	Oxygen Carrying Bag WaterStop PROFI (PLANE)	approx. 16,4 L	approx. 49 x 24 x 17 cm	80,00 €**



## Ordering information
























Image	Unit	PU	Art. No.	Description	Tragkraft/ Volume	Measures	SRP
	Piece	10	N190312 (red)	Rescue sheet PRO SAFE (PLANE)	bis 150 kg	200 x 70 cm	55,00 €**
	Piece	10	N190206 (red)	Rescue sheet PRO (PLANE)	bis 150 kg	200 x 70 cm	38,00 €**
	Piece	10	N190457 (Black)	Rescue sheet PRO (PLANE)	bis 150 kg	200 x 70 cm	40,00 €**
<b>NEW</b> 	Piece	10	N190458 (Black)	Rescue sheet <b>ONE WAY</b> BLACK L	bis 200 kg	200 x 80 cm	15,00 €**
<b>NEW</b> 	Piece	10	N190459 (Black)	Rescue sheet <b>ONE WAY</b> BLACK XL	bis 350 kg	200 x 98 cm	20,00 €**
	Piece	4	N180035 (red) N180064 (green on request) N180065 (yellow)	Emergency Bag FLEX (PLANE)	approx. 17,3 L	approx. 28 x 32 x 18 cm	60,00 €**
	Paar	50	N180044	Wallhooks for N180035, N180064, N180065	--	approx. 7 x 2,5 x 4,5 cm	6,00 €**
	Piece	20	N180043 (red) N180042 (Blue) N180041 (red) N180040 (Blue)	Driver's logbook A5 Organizer (PLANE)	approx. 3,1 L	approx. 28 x 5 x 22 cm	25,00 €**
	Piece	20	N180038 (red) N180037 (Blue)	Operation Organizer A4 MULTI (PLANE)	approx. 6,2 L	approx. 28 x 6 x 37 cm	48,00 €**
	Piece	2	N190189-V2 (red) N190171-V2 (Blue) N190441 (green on request) N190442 (yellow)	PREMIUM Emergency Backpack X1 YELLOW (PLANE)	approx. 39 L	approx. 52 x 30 x 25 cm	200,00 €**
	Piece	2	N190115 (red) N190117 (Blue)	Emergency Backpack PROFI (PLANE)	approx. 67,1 L	approx. 65 x 42 x 23 cm	160,00 €**
	Piece	2	N190114 (red)	Emergency Backpack PROFI (NYLON)	approx. 67,1 L	approx. 65 x 42 x 23 cm	120,00 €**
	Piece	4	N190134 (red) N190244 (Blue)	Emergency Backpack YELLOW LARGE (PLANE)	approx. 28 L	approx. 36 x 47 x 26 cm	69,00 €**
	Piece	10	N190132 (red) N190243 (Blue) N190133 (Black) N190445 (green on request) N190446 (yellow)	Emergency Backpack YELLOW/ WHITE LARGE (NYLON)	approx. 28 L	approx. 36 x 47 x 26 cm	48,00 €**
	Piece	6	N190137 (red) N190241 (Blue)	Emergency Bag YELLOW LARGE (PLANE)	approx. 28 L	approx. 52 x 32 x 32 cm	69,00 €**
	Piece	10	N190135 (red) N190242 (Blue) N190136 (Black)	Emergency Bag YELLOW LARGE (NYLON)	approx. 28 L	approx. 52 x 32 x 32 cm	48,00 €**
	Piece	10	N190139 (red) N190239 (Blue)	Emergency Bag YELLOW MEDIUM (PLANE)	approx. 18 L	approx. 43 x 25 x 27 cm	45,00 €**
	Piece	20	N190138 (red) N190240 (Blue) N190435 (green on request) N190436 (yellow)	Emergency Bag YELLOW MEDIUM (NYLON)	approx. 18 L	approx. 43 x 25 x 27 cm	25,00 €**
	Piece	12	N190140 (red) N190203 (Blue) N190447 (green on request) N190448 (yellow)	Emergency Bag YELLOW SMALL (NYLON)	approx. 12 L	approx. 37 x 21 x 27 cm	25,00 €**

Image	Unit	PU	Art. No.	Description	Volume	Measures	SRP
	Piece	100	N190370 (Blue) N190371 (green) N190372 (red) N190373 (yellow) N190374 (Black / N190375 (grau)	WaterStop Module bag PROFI XS (PLANE)	approx. 1,4 L	approx. 16 x 15 x 6 cm	7,50 €**
	Piece	100	N190377 (Blue) N190378 (green) N190379 (red) N190380 (yellow) N190381 (Black) N190382 (grau)	WaterStop Module bag PROFI S (PLANE)	approx. 2,3 L	approx. 22 x 15 x 7 cm	7,80 €**
	Piece	100	N190384 (Blue) N190385 (green) N190386 (red) N190387 (yellow) N190388 (Black) N190389 (grau)	WaterStop Module bag PROFI M (PLANE)	approx. 1,8 L	approx. 25 x 12 x 6 cm	8,00 €**
	Piece	50	N190391 (Blue) N190392 (green) N190393 (red /PLANE) N190394 (yellow) N190395 (Black) N190396 (grau)	WaterStop Module bag PROFI L (PLANE)	approx. 5,2 L	approx. 37 x 14 x 10 cm	8,70 €**

## Suitable for:







HERZmed WaterStop  
Emergency Backpack  
PROFI



HERZmed WaterStop  
Emergency Backpack  
PRO



HERZmed Emergency  
Backpack  
YELLOW Large

Image	Unit	PU	Art. No.	Description	Volume	Measures	SRP
	Piece	50	N180058	Intubationstasche MULTI (NYLON)		approx. 33 x 18 x 3 cm	25,00**
	Piece	20	N180059	Intubationstasche ULTRA (PLANE)	for approx. 100 Ampullen	approx. 29 x 24 x 8,5 cm	35,00**
	Piece	20	N180060 (Blue)	Aufbewahrungstasche Videolaryngoskop blue (PLANE)		approx. 29 x 24 x 8,5 cm	45,00**
	Piece	50	N190077 (red)	Intubationstasche RED (PLANE)		approx. 34 x 23 x 6 cm	15,00**

## NEW IN THE ASSORTMENT!

	Piece	4	N180200 (red) N180201 (Blue) N180202 (Black) N180203 (yellow)	Emergency Bag AED FLEX (PLANE)	approx. 30 L	approx. 44 x 30 x 25 cm	95,00**
	Piece	2	N190434 (yellow)	HERZmed Emergency Bag Yellow L PLANE	approx. 32 L	approx. 50 x 30 x 23 cm	69,00**



## Ordering information

Image	Unit	PU	Art. No.	Description	Volume	Measures	SRP
	Piece	50	N190124 (Blue)	Module bag S (PLANE)	approx. 2,6 L	approx. 24 x 12 x 9 cm	5,50 €**
	Piece	50	N190127 (Blue) N190128 (red) N190129 (yellow) N190130 (orange) N190131 (green on request)	Module bag S (NYLON)	approx. 2,6 L	approx. 24 x 12 x 9 cm	4,50 €**
	Piece	50	N190123 (red, PLANE) N190125 (green on request, NYLON)	Module bag M	approx. 1,9 L	approx. 32 x 12 x 5 cm	5,50 €**
	Piece	50	N190122 (green on request)	Module bag L (PLANE)	approx. 4,6 L	approx. 23 x 20 x 10 cm	5,50 €**
	Piece	50	N190126 (Blue)	Module bag L (NYLON)	approx. 4,6 L	approx. 23 x 20 x 10 cm	4,80 €**
	Piece	100	N180107 (Blue) N180108 (green on request) N180109 (red) N180110 (yellow) N180111 (Black)	WaterStop Module bag S (PLANE)	approx. 2,6 L	approx. 24 x 12 x 9 cm	6,50 €**
	Piece	100	N180086 (Blue) N180087 (green on request) N180088 (red) N180089 (yellow) N180090 (Black)	WaterStop Module bag S (NYLON)	approx. 2,6 L	approx. 24 x 12 x 9 cm	6,00 €**
	Piece	100	N180093 (Blue) N180094 (green on request) N180095 (red) N180096 (yellow) N180097 (Black)	WaterStop Module bag M (NYLON)	approx. 1,9 L	approx. 32 x 12 x 5 cm	6,60 €**
	Piece	100	N180100 (Blue) N180101 (green) N180102 (red) N180103 (yellow) N180104 (Black)	WaterStop Module bag L (NYLON)	approx. 4,6 L	approx. 23 x 20 x 10 cm	6,80 €**

### Suitable for:



HERZmed WaterStop  
Emergency Backpack  
PROFI



HERZmed WaterStop  
Emergency Backpack PRO



HERZmed Emergency  
Backpack  
YELLOW Large



HERZmed WaterStop  
Emergency Backpack  
SMALL



HERZmed WaterStop  
Emergency Bag FLEX

Image	Unit	PU	Art. No.	Description	Volume	Measures	SRP
	Piece	50	N190349 (Blue) N190350 (green on request) N190351 (red) N190352 (yellow) N190354 (grau)	WaterStop Module bag X1 S (PLANE)	approx. 3,2 L	approx. 21 x 15 x 10 cm	11,00 €**
	Piece	50	N190356 (Blue) N190357 (green on request) N190358 (red) N190359 (yellow) N190361 (grau)	WaterStop Module bag X1 L (PLANE)	approx. 3,9 L	approx. 30 x 13 x 10 cm	10,00 €**

### Suitable for:



HERZmed PREMIUM  
Emergency Backpack X1

**NEW IN THE ASSORTMENT!**

## HERZmed® Emergency Bag AED FLEX PLANE

## HERZmed® Emergency Bag Yellow L PLANE



### Technical specifications

The technical specifications can be found with the respective products. If you have any questions, please contact our customer service.

**The material makes the difference** WaterStop® – for demanding applications

We are pleased to introduce our new WaterStop product line. We have been supplying a very successful range of emergency containers for years. In addition, we are keen to improve our existing range through innovation, higher quality high-tech materials and manufacturing processes. Our aim is to respond to the existing requirements and needs of professional users.

**WaterStop products are designed for professional use.**

#### Water-repellent zip fasteners

The particularly robust and special zips are water-repellent and therefore also suitable for use in the rain and enable the storage of dressing material that is susceptible to moisture.

- Water-repellent material
- Easier cleaning of conventional soiling
- Improved service life
- Easy opening and closing



#### PU-coated material

The polyester fabric is coated twice with a high-quality PU coating. This has numerous advantages over the conventional PVC coating.


- Significantly lower weight
- Higher load capacity in daily use
- Better temperature resistance, especially at low temperatures the material remains flexible
- Improved service life
- More environmentally friendly
- than conventional PVC material



#### WaterStop® Technology with TEFLON®

Polyester fabric with TEFLON® coating, which turns ordinary polyester into a unique material. Water and oil can no longer penetrate the fabric and are repelled on the surface.

- Water-repellent material
- Improved durability
- Easier cleaning of conventional soiling
- Significantly reduced weight



### Available materials

#### Available colours



Would you like to customise one of our models individual requirements? Our service staff are at your disposal.

Tel.: 0911 - 14 88 959 - 0  
Mail: fachhandel@bexamed.de

### Ordering and counselling by:

HERZmed GmbH

Herkulan-Schwaiger-Gasse 16  
82487 Oberammergau

Tel.: +49 (0) 911 / 14 88 959 - 0  
Fax: +49 (0) 911 / 14 88 959 - 9

Mail: herzmed@herzmed.de  
Web: www.herzmed.de

The products shown are only an excerpt from our product range. Exact details and product information are available on request or at [www.herzmed.de](http://www.herzmed.de). All prices marked with an asterisk are plus statutory VAT. All brand names used remain the property of the respective brand owner. For errors, possible printing errors and colour deviations. Our General Terms and conditions apply, which you can also view at [www.herzmed.de](http://www.herzmed.de).



Scan for German version  
[herzmed.com/downloads](http://herzmed.com/downloads)

