

Rescue sheet

ONE-WAY BLACK



Carries more than responsibility

The HERZmed ONE WAY rescue sling is suitable for transporting people from rough terrain or areas that are difficult to access. It can be used, for example, by rescue services, the fire department or mountain rescue services.

Load straps, which are sewn crosswise to the sling, provide high stability. The HERZmed ONE WAY rescue sling is used for both seated and prone transportation of patients, whereby they are transported as if wrapped in a carrier bag.

The HERZmed ONE WAY rescue sling is made of high-strength polyester fabric coated on both sides with PVC. The load straps are made of woven polyester covered with PE plastic tubing.

Details



Additional protection with headboard



and footmuff



double PVC-coated polyester tarpaulin material

Highlights

- With leg bag and head protection
- Water-resistant and hard-wearing PVC coated polyester tarpaulin material
- 10/14 handles
- Size of the lying surface: L 200 cm x 80 cm
- Size of the lying surface: XL 200 cm x 98 cm
- Load capacity L 200 kg or XL 250 kg
- Rescue sheet complies with DIN EN 1865
- Can also be used as a rescue seat

Lightweight and compact

Space-saving and lightweight.

Headboard

Integrated headboard for protection.

High load capacity

Load capacity of up to
L = 200 kg
XL = 250 kg

Size L 10 / XL14 Straps

The straps are made of woven polyester material covered with PE plastic tubing.

Large lying surface

Version L = 200 x 80 cm
Version XL = 200 x 98 cm.

Durable

Made from hard-wearing coated on both sides with PVC coated high-strength polyester fabric.

Impermeable to liquids and 100% X-ray transparent.

Leg bag

The leg bag prevents the patient from slipping out.



Versatile in use

Load straps, which are sewn crosswise to the sling, provide high stability. The HERZmed One Way rescue sling can be used for transporting patients in a sitting or lying position, whereby they are transported as if wrapped in a carry bag. Disposal: In normal household waste as residual waste.

Available in:

Size L



Size XL



Technical specifications

Size:

Version L = 200 x 80 cm
Version XL = 200 x 98 cm

Load capacity:

Version L = 200 kg
Version XL = 250 kg

Colour: BLACK

Packaging:

Individually packed in transparent foil

Sterile: No

Storage: Store in a dry place

Disposal: In normal household waste as residual waste

Materials:

Lying surface:
Made of hard-wearing high-strength polyester fabric coated on both sides with PVC. Impermeable to liquids and 100% X-ray transparent.

Carrying straps made of durable woven polyester covered with PE plastic tubing.

Order information:

| Art.-No. | Description | Image |
|----------|--------------------------------|---|
| 190458 | Rescue sheets ONE-WAY BLACK L |  |
| 190459 | Rescue sheets ONE-WAY BLACK XL |  |

Order and advice:

Vertrieb/Sales Office:
HERZmed GmbH
Rottenbacher Straße 9
82497 Unterammergau

Tel.: +49 (0) 8822 824 0002
Fax: +49 (0) 8822 824 0006

Verwaltung/Management:
HERZmed GmbH
Herkulan-Schwaiger-Gasse 16
82487 Oberammergau

Mail: herzmed@herzmed.de
Web: www.herzmed.com

The products shown are only an excerpt of our product range. Exact details and product information are available on request or at www.herzmed.de. All prices marked with an asterisk are exclusive of VAT. All brand names used remain the property of the respective brand owner. We assume no liability for errors, possible misprints and color deviations. Our general terms and conditions apply, which you can also view at www.herzmed.de.

